



# PIPRO

Finger  
Extension  
Orthosis

## Wearing Your PIPRO

Great outcomes start with a properly fitted splint from your therapist. Your wear schedule:

1. \_\_\_\_\_ minutes wearing \_\_\_times/day for \_\_\_\_ days.
2. If your splint is comfortable and skin looks fine after wearing the splint, progress to \_\_\_\_\_ minutes wearing during the day for several days.
3. If comfortable, skin in feels and looks good and your therapist recommends. progress to night time wearing.



## Cleaning

Hand wash with mild soap/detergent in cool water. Rinse thoroughly. Air dry only. Do not bleach.



## Precautions

*This Product Contains Natural Rubber Latex which may cause allergic reactions.*

Stop using your splint and call your health professional if you develop skin redness, a blister, numbness, coolness or pain in your finger



Visit [angularortho.com](http://angularortho.com) for further information and videos. It's easy by clicking the OR code